

Helps the patient to have a bath or shower without getting their PICC line, midline or elbow dressing wet.



Prescription Request Form

Patient _____

Hospital / Clinic _____

Clinician _____

Date _____

Elbow LimbO – Seals above and below the elbow, ideal protection for PICC lines, midlines or elbow dressings.
Sizing instructions: Check your patient's weight and arm circumference to select the correct size elbow LimbO.



Patient Weight	Arm Circumference above PICC Line or Dressing	Model	Size	PIP Code	✓
51 to 70kg (8 – 11 Stone)	25 - 29cm (10 - 11.5")	MP65	Slim / Medium	386 - 4303	
70kg + (11 Stone +)	30 - 39cm (12 - 15")	MP75	Medium / Large	386 - 4311	

Advice line: **01243 578755** Email: info@limboproducts.co.uk www.limboproducts.co.uk

Thesis Technology Products Ltd., Brooks Green Farm, Brooks Lane, Bosham, West Sussex, PO18 8JX

EPP0222