Helps the patient to have a bath or shower without getting their PICC line, dressing or compression hosiery wet.



Prescription Request Form

T delotte	1100011017		
Clinician) 0:00	Date	
Half Leg LimbO - Seals just above	the knee, ideal for lower leg bar	ndages, casts or bulky foot dressings	. For bath or shower.

Hospital / Clinic

Sizing Instructions: Check your patient's height and weight to select the correct size half leg LimbO.



Patient

Patient Height	Patient Weight	Model	Description	PIP Code	\checkmark
5'5" (165 cm)	Up to 10 Stone (64kg)	MP76	Slim Leg	333-7581	
and above	10 to 16 Stone (64 to 102kg)	MP80	Average Build	320-0045	
	16 to 22 Stone (102 to 140kg)	MP180	Large Leg	320-0060	
Under 5'5" (165cm)	Up to 10 Stone (64kg)	MP76S	Slim Leg Short	333-7599	
(======,	10 to 16 Stone (64 to 102kg)	MP80S	Short Leg	320-0052	
	16 to 22 Stone (102 to 140kg)	MP180S	Large Leg Short	320-0078	
	140Kg+ (22 Stone+)	MP200HL	Extra Large Leg	431-9224	

Foot LimbO - Seals just above the ankle, ideal for light dressings on the foot. Shower use only.



Leg Circumference (5cm above ankle bone)	Model	Size	PIP Code	\checkmark
20 - 25cm (8-10")	MP20	Small	416-6476	
23 - 34cm (9-13")	MP25	Medium to Large	416-6484	

Please note: Maximum UK shoe size is 11. If larger, consider a half leg LimbO option.

Advice line: 01243 573417 Email:info@limboproducts.co.uk www.limboproducts.co.uk

PLEASE TURN OVER FOR ELBOW, ARM AND FULL LEG OPTIONS.