Helps the patient to have a bath or shower without getting their PICC line, dressing or compression hosiery wet.



## **Prescription Request Form**

Patient			_ Hospital / Clinic		
Clinician		0.00		Date	
70			and the second		

Half Leg LimbO - Seals just above the knee, ideal for lower leg bandages, casts or bulky foot dressings. For bath or shower. Sizing Instructions: Check your patient's height and weight to select the correct size half leg LimbO.



	Patient Height	Patient Weight	Thigh	Model	Description	PIP Code	$\checkmark$
	5'5" (165 cm) and above	Up to 10 Stone (64kg)	35-39cm	MP76	Slim Leg	333-7581	
		10 to 16 Stone (64 to 102kg)	39-54cm	MP80	Average Build	320-0045	
		16 to 22 Stone (102 to 140kg)	52-57cm	MP180	Large Leg	320-0060	
	Under 5'5" (165cm)	Up to 10 Stone (64kg)	35-39cm	MP76S	Slim Leg Short	333-7599	
		10 to 16 Stone (64 to 102kg)	39-54cm	MP80S	Short Leg	320-0052	
		16 to 22 Stone (102 to 140kg)	52-57cm	MP180S	Large Leg Short	320-0078	
		140Kg+ (22 Stone+)	65-90cm	MP200HL	Extra Large Leg	431-9224	

Foot LimbO - Seals just above the ankle, ideal for light dressings on the foot. Shower use only.



Leg Circumference (5cm above ankle bone)	Model	Size	PIP Code	$\checkmark$
20 - 25cm (8-10")	MP20	Small	416-6476	
23 - 34cm (9-13")	MP25	Medium to Large	416-6484	

Please note: Maximum UK shoe size is 11. If larger, consider a half leg LimbO option.

Advice line: 01243 573417 Email:info@limboproducts.co.uk

www.limboproducts.co.uk

PLEASE TURN OVER FOR ELBOW, ARM AND FULL LEG OPTIONS.